

# Zeronil 402 mg Spot-on Solution for extra large dogs

Authorised

- Fipronil

## Product identification

**Medicine name:**

Zeronil 402 mg Spot-on Solution for extra large dogs

---

**Active substance:**

Fipronil

---

**Target species:**

Dog

---

**Route of administration:**

Spot-on use

---

## Product details

**Active substance and strength:**

Fipronil

402.00 milligram(s) / 1.00 Pipette

---

**Pharmaceutical form:**

Spot-on solution

---

**Anatomical therapeutic chemical veterinary (ATCvet) codes:**

QP53AX15

---

**Legal status of supply:**

Veterinary medicinal product not subject to veterinary prescription

---

**Authorisation status:**

Valid

---

**Authorised in:**

Hungary

---

**Package description:**

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 150 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 90 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 60 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 30 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 24 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 21 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 20 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 18 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 15 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 12 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 10 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 9 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 8 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 6 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer.A Cardboard box with 4 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer. A Cardboard box with 3 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer. A Cardboard box with 2 pipettes in individual foil sachets.

A 4.02ml white pipette composed of a heat-formed shell of a polypropylene/cyclicolefin copolymer/polypropylene layer and a polyethylene/ethylene vinylalcohol/polyethylene layer. A Cardboard box with 1 pipette in individual foil sachets.

---

## Additional information

**Entitlement type:**

Marketing Authorisation

---

**Legal basis of product authorisation:**

Hybrid application (Article 13(3) of Directive No 2001/82/EC)

---

**Marketing authorisation holder:**

Chanelle Pharmaceuticals Manufacturing Limited

---

**Marketing authorisation date:**

26/03/2012

---

**Manufacturing sites for batch release:**

Chanelle Pharmaceuticals Manufacturing Limited

---

**Responsible authority:**

Directorate Of Veterinary Medicinal Products

---

**Authorisation number:**

3122/X/12 NÉBIH ÁTI

---

**Date of authorisation status change:**

26/03/2012

---

**Reference member state:**

Ireland

---

**Procedure number:**

IE/V/0276/004

---

**Concerned member states:**

Bulgaria Czechia Estonia France Germany Greece Hungary Italy Latvia  
Liechtenstein Lithuania Luxembourg Malta Netherlands Poland Portugal  
Romania Slovakia Spain Sweden

---

To consult adverse reactions on veterinary medicinal products please go to  
[www.adrreports.eu/vet](http://www.adrreports.eu/vet)

## Documents

Summary of Product Characteristics

ie-puar-mr-iev0276004-zeronil-402-mg-spot-on-solution-for-extra-large-do-en.pdf